
SVBONY

EN

DE

FR

IT

ES

RU

JP

SM201 Series
Biological Microscope
General Instructions

Contents

Applications of Microscope	-----01
Working Principle	-----01
Structure & Main Technical Specification	-----01
Installation and Operation	-----03
Microscope Maintenance	-----03
Guarantee	-----29

EN(01-04)

DE(05-08)

FR(09-12)

IT(13-16)

ES(17-20)

RU(21-24)

JP(25-28)

EU Importer

Name: Germany Retevis Technology GmbH

Address: Uetzenacker 29,38176 wendeburg

1. Applications of Microscope:

This series of microscopes can be used for biological, pathological, bacteriological observation, teaching, professional research, clinical experiment and routine medical examination in medical and health institutions, laboratories, agricultural science institutes, research institutes and colleges and universities. This product is easy to use, safe and reliable. With appropriate auxiliary equipment, the microscope can play a huge advantage in the demonstration of education, image acquisition, storage and analysis of processing.

2. Working Principle

The working principle of the microscope is shown in Figure 1. The illuminator (a) provides the light source, and the light passes through the condenser lens system (b) and becomes suitable light to illuminate the specimen sheet (c), after changing the direction of transmission, it passes through the eyepiece system (f) and zooms in for the second time into the viewer's field of vision.

The total magnification of the microscope = (the magnification of the objective lens in use) X (the magnification of the eyepiece in use)


Figure 1. Schematic diagram of microscope

3. Structure & Main Technical Specification

The structure is as shown in Figure 2.

1. Eyepiece

The eyepieces of this microscope are routinely configured as wide-field eyepieces WF10X and WF20X. Figure 1. Schematic diagram of microscope

The basic specifications are as follows:


Wide-field eyepieces	Magnification	Field Diameter	Working distance	Remark
WF10X	10	15mm	24.95mm	
WF25X	25	8mm	17mm	

2. Objective Lens

The standard configuration of the objective lens of this microscope is a DIN achromatic objective lens, mainly including 4X, 10X, 40X (spring), 100X (spring head, Oil immersion lens). The 100X objective lens is an oil immersion lens. When using it, immersion oil should be added between its front lens and the cover glass. If there are air bubbles entering the oil layer, the image quality will deteriorate. To drive out the air bubbles, you can turn the converter several times or add a little more oil. After the oil immersion observation, remember to clean the front lens of the objective lens and the parts stained with oil one by one.

The basic parameters of the DIN achromatic objective lens system are as follows:

DIN Achromatic Objective Lens	Magnification	Aperture	Focal Length	Working Distance	Coverslip Thickness	Remark
4X	4	0.10mm	31.04mm	37.5mm	0.17mm	
10X	10	0.25mm	17.13mm	7.316mm	0.17mm	
40X	40	0.65mm	4.65mm	0.632mm	0.17mm	Spring
100X	100	1.25mm	2.906mm	0.198mm	0.17mm	Spring head, Oil immersion lens

- | | | | |
|--------------------------|---------------------|-------------|-------------------------|
| 1. Eyepiece | 2. Observation head | 3. Switch | 4. Objective lens |
| 5. Platform | 6. condenser | 7. Rubber | 8. Cable |
| 9. Base | 10. Arm | 11. Focuser | 12. Platform limit nuts |
| 13. Platform limit screw | | | |

Figure2 Microscope Structure

3. Mechanical Strength: 160mm

4. Objective conjugate distance: 195mm

5. Observation head:

Observation head is mainly used to change the direction and path of light propagation to adapt to different users and achieve a variety of observation effects. This microscope is standard equipped with binocular observation head, 45° tilt, 360° rotation.

6. Converter:

The converter adopts advanced rolling bearing structure, which is well-processed and positioned accurately, which can better meet the requirements of microscope alignment and center unity, making the operation comfortable and convenient for the observer.

7. Platform:

Double-layer mechanical mobile platform, with N.A.1.25 Abbe condenser, color filter holder, rainbow iris diaphragm.

8. Focus System:

The microscope adopts a coaxial rack and pinion focusing mechanism with coarse and fine movements. The focus range of coarse and fine movement is 14mm, and the fine adjustment step value is 0.015mm. Turn the coarse focusing handwheel (large), the platform (5) will rise and fall quickly; turn the fine handwheel (small), the platform (5) will slowly and accurately rise and fall. The platform limit nut (12) and the platform limit screw (13) limit the rising limit position of the platform (5), so that the lifting of the platform (5) has absolute safety and reliability during focusing, and will not be caused by the platform (5) The rise is too high and the objective lens breaks the glass section. In general, it has been adjusted when the product leaves the factory. In case of deviation during use, it can be adjusted as follows: first loosen the platform limit nut (12), screw up the platform limit screw (13), put the 40X objective lens at the working position, and then place it on the platform (5) Put the specimen slice on top, then raise the platform (5) to the position where the specimen slice contacts the objective lens, keep the platform still, turn down the platform limit screw (13) to contact the guide rail in the rack, and use the platform limit nut (12) to secure the position of the platform stop screw (13).

9. Light Source System

LED light, soft light, long life (about 100,000 hours), good safety. It generates less heat and won't burn your hands no matter how long you use it. Equipped with a rechargeable battery, it can be used in a working environment without electricity temporarily. The center of the illuminator is important to the imaging quality of the microscope. When the center of the illuminator deviates from the center of the microscope objective lens, there may be a bright side and a dark side in the field of view. The illuminator has a certain service life. When replacing the illuminator, it must be consistent with the illuminator equipped with the instrument. Before changing the illuminator, unplug the power plug of the instrument. After the microscope cools down, put the microscope backwards, take out the illuminator door cover screw (15), flip up the door cover, take out the old illuminator, and replace it with a new illuminator to ensure that the illuminator is safe. The center is in the correct position. Return the door cover to its original position, and tighten the illuminator door cover screws (15).

Do not touch the glass part of the illuminator with your hands. Wear gloves or cover the illuminator with a protective cover when installing. Use a clean cloth dipped in pure alcohol to wipe off fingerprints and smudges on the glass, because these things on the surface of the illuminator will weaken it. brightness and cause the illuminator to burst.

During lighting or just after lighting, the illuminator and heat dissipation ribs will be very hot. Be careful not to let the illuminator burn yourself when changing the illuminator, but don't panic.

4. Installation and Operation

1. Open the package carefully, and count all the parts in it according to the packing list.

2. For convenience of packing and shipping, some of the components are unassembled with the main product, and please assemble them according to Figure 2 before use.

3. Operation steps:

(1) Select a working environment with no direct light on the device, keep the device away from a window and don't face a window directly, because the direct light will affect the contrast of images and observation results. Working environment of the telescope:

- a) Indoor temperature: 0 C~40 C; Maximum Relative Humidity: 85%
- b) High temperature and relative humidity will cause mildew, condensation and damage to microscope.
- c) Avoid placing the microscope in a dusty environment. When the microscope is not in use, it should be covered with a microscope dust cover.
- d) The microscope should be stably placed where it is flat and with no vibration.

(2) Carefully clamp the specimen slice in the clamps. When loading the slices, ensure that the cover glass is facing the objective lens, otherwise, no image can be observed no matter how the focus is adjusted when using a high-magnification objective lens. Turn the low-magnification (4X or 10X) to the working position, lift the platform until it touches the objective lens, and check whether the limited lifting position is correct.

(3) Turn the low-power objective (4X or 10X) to the working position and turn on the power (or adjust the reflector Angle).

Attention: The supply voltage should be consistent with the calibration voltage of the microscope, otherwise it will lead to line damage and safety accidents.

(4) Move the specimen slice so that the center of the slice is at the center of the condenser.

(5) Use the coarse handwheel (large) to adjust the focal length of the objective lens until the image of the specimen is bright and clear. When operating, generally use a low-magnification objective lens to determine the approximate position of the focal plane, and use this to understand the focusing direction of the focusing handwheel. Finally, use the micro handwheel (small) to fine-tune the image to obtain a satisfactory clear image.

(6) After the fine focus is completed, if you turn the converter (3) to change other objective lenses for observation, you can see the image at this time, but it may not be very clear, just adjust the fine-tuning handwheel (small) a little bit, this is because All objectives on a microscope are parfocal and maintain the same center.

Notice:

Although the 4X and 10X objective lenses are integrally fixed, because of their long focal lengths, they will not break the specimen slices during observation.

The focal lengths of 40X and 100X are very short, and the front lens is almost attached to the specimen cover glass during observation, but because they have a spring device inside, they can expand and contract, so they will not easily break the specimen slice.

(7) In order to observe a clearer image, different diameters of the variable diaphragm (dial diaphragm or rainbow diaphragm) can be selected according to the objective lens used.

5. Microscope Maintenance

1. Be careful when opening the box to prevent the lens and other accessories from falling and being damaged.

2. All lenses have been assembled and adjusted, please do not disassemble them by yourself.

3. The structure of the objective lens converter and the coarse and fine focus adjustment mechanism is compact, please do not disassemble and assemble easily, if the instrument breaks down, please send it to a professional repair shop for repair.

4. The instrument should be kept clean, and the dust should be removed frequently. The sliding parts should be regularly coated with a small amount of non-corrosive grease. When cleaning, special attention should be paid not to pollute the optical parts.

5. The instrument should be placed in a cool and dry place. After use, the power must be cut off and cover the microscope with a microscope dust cover. If it is not used for a long time, it is best to unscrew the objective lens and put it in the mirror box, and screw the dust-proof cover on the converter.

Before using this device, read this guide which contains important operating instructions for safe usage and control for compliance with applicable standards and regulations.

FCC Requirements:

- Products authorized under Part 15 using SDoC or Certification require a label containing one of the following compliance statements

- (1) Receivers associated with licensed device service operations:
This device complies with part 15 of the FCC Rules. Operation is subject to the condition that this device does not cause harmful interference.
- (2) Stand-alone cable input selector switch:
This device complies with part 15 of the FCC Rules for use with cable television service.
- (3) All other devices:
 - This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
 - (1) This device may not cause harmful interference, and
 - (2) this device must accept any interference received, including interference that may cause undesired operation.

CE Requirements:

- (Simple EU declaration of conformity) Hong Kong Svbon Technology Co., Ltd declares that the equipment type is in compliance with the essential requirements and other relevant provisions of RED Directive 2014/30/EU and the ROHS Directive 2011/65/EU and the WEEE Directive 2012/19/EU; the full text of the EU declaration of conformity is available at the following internet address: www.svbon.com.

• Disposal

The crossed-out wheeled-bin symbol on your product, literature, or packaging reminds you that in the European Union, all electrical and electronic products, batteries, and accumulators (rechargeable batteries) must be taken to designated collection locations at the end of their working life. Do not dispose of these products as unsorted municipal waste. Dispose of them according to the laws in your area.

IC Requirements:

CAN ICES-3(B)/NMB-3(B)

Avoid Choking Hazard

Small Parts. Not for children under 3 years.

Approved Accessories

WARNING

- This device meets the regulatory standards when used with the Svbon accessories supplied or designated for the product.
- For a list of Svbon-approved accessories for your item, visit the following website: <http://www.Svbon.com>.

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Post Code: _____ Email: _____

Remarks:

1. This guarantee card should be kept by the user, no replacement if lost.
2. Most new products carry a one-year manufacturer's warranty from the date of purchase.
3. The user can get warranty and after-sales service as below:
 - Contact the seller where you buy.
4. For warranty service, you will need to provide a receipt proof of purchase from the actual seller for verification

Exclusions from Warranty Coverage:

1. To any product damaged by accident.
2. In the event of misuse or abuse of the product or as a result of unauthorized alterations or repairs.
3. If the serial number has been altered, defaced, or removed.

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Hong Kong Svbon Technology Co.,Ltd

Add: Unit B, 5th Floor, Gallo Commercial Building, 114-118
Lockhart Road, Wanchai, Hong Kong

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