
SVBONY

SM402
Zoom-stereo microscope
User Manual

CONTENTS

1.Matters need attention.....	1
2.Name of each component.....	1
3.Installation of various parts.....	1
4.Operation.....	2
5.Technical Parameters.....	2
6.Fault handling.....	3

DE (5-8)

FR (9-12)

IT (13-16)

ES (17-20)

JP (21-24)

EU Importer: Germany Retevis Technology GmbH
Address: Uetzenacker 29,38176 wendeburg

1. Matters need attention


Operation

- (1) The components should not be directly exposed to the sun, and should be placed in a dry and clean environment to avoid high temperature and severe vibration.
- (2) The microscope is a precision instrument. It should be handled with care and avoid impact and collision during transportation.
- (3) In order not to affect the clarity of the image, avoid dirt or fingerprints on the surface of the lens.
- (4) Need to hold the focusing wheel with both hands, rotate and adjust in the same direction, otherwise it will malfunction.


Care and Maintenance

- (1) All lenses should be kept clean. If there is fine dust, it can be blown off with a balloon or gently wiped off with cotton yarn; if there are oil stains and fingerprints, it can be wiped off lightly with cotton yarn dipped in a small amount of ethanol and ether mixture at a ratio of 3:7.
- (2) Do not use organic solvents to wipe the remaining surfaces of the microscope, especially the surface of plastic products, and neutral detergent should be used for cleaning.
- (3) Do not disassemble or assemble the microscope by yourself, so as not to affect the performance of the microscope.
- (4) When the microscope is not in use, it should be covered with a dust cover so that it does not come into contact with dust and stored in a place isolated from moisture to avoid rust or mold.
- (5) In order to maintain the performance of the microscope, it is recommended to conduct regular inspections (for details, please contact the nearest agent).

2. Name of each component


3. Installation of various parts

4. Operation

4-1 Use of glass work board

(1) You can lift the work board directly from the recess with your fingertips (Figure 1).


(Figure 1)

4-2 Adjustment of the tightness of the focusing mechanism

(1) To adjust the tightness of the focusing mechanism, you can hold one of the handwheels with your hand and rotate the other handwheel. The tightness depends on the direction of rotation of the handwheel; tighten it clockwise and loosen it counterclockwise (see Figure 2).

(2) Adjusting the tightness of the focusing mechanism properly can prevent the microscope from sliding down with the bracket during the observation of the microscope body, and also make Focusing more comfortable.

(Figure 2)

4-3 Specimen placement

(1) Place the specimen in the middle of the glass work plate, and if necessary, press the specimen with a press sheet.

(2) Use the illuminator to illuminate the specimen.

4-4 Diopter adjustment and focus


(1) Rotate the zoom handwheel to the maximum magnification.

(2) Rotate the diopter adjustment ring to the 0 mark position.

(3) Observe through the eyepiece on the right. If the image is not clear, turn the focusing handwheel to make the specimen image clear.

(4) Rotate the zoom handwheel to the minimum magnification.


(5) Observe through the eyepiece on the right, if the image is not clear, rotate the right diopter adjustment ring ② to make the specimen clear (Figure 3)

(Figure 3)

(6) Then rotate the zoom handwheel to the maximum magnification. Observe through the eyepiece on the right, if the image is not clear, repeat the above 3 to 5 steps, which can make the diopter adjustment more precise.

(7) Rotate the zoom handwheel to the minimum magnification. Observe through the left eyepiece, if the image is not clear, rotate the left diopter adjustment ring ① to make the image clear (Figure 3).

(Figure 4)

4-5 Adjustment of interpupillary distance


(8) Hold the left and right prism boxes with your hands, and turn the left and right prism boxes in the direction shown by the arrows in (Figure 4)

until the binocular observation is comfortable.

4-6 Use of eyepiece cover

(1) For people who do not wear glasses, hold the diopter adjustment ring with your hands so that it does not rotate. By rotating the eyepiece, the eyepiece cover can fit the observer's eye better.

(2) For people who wear glasses, the eyepiece cover can be removed for observation.

(Figure 5)


4-7 Installing and removing the reticle

(1) Unscrew the pressure ring from the eyepiece (as shown in the figure).

(2) Clean the reticle ①. Place the engraved side of the reticle facing down and fit it into the press ring (Figure 5).

(3) Screw the pressure ring ② with the reticle into the eyepiece until it is tightened.


(4) To remove the reticle, unscrew the pressing ring ③ from the eyepiece, Take out the reticle, wrap it in clean paper, and store it.

(Figure 6)

4-8 Installation of luminaire

Insert the boss of the illuminator assembly into the ring LED bracket facing the objective lens, and then tighten the screws (as shown in Figure 6).


(Figure 7)

4-9 Choose light path

(1) Use the lever to move inward and outward to switch between binocular observation and photography. Move the lever outwards for photography and video (When selecting photography and recording, visual inspection can also be carried out at the same time); conversely, move the lever inwards for binocular observation.

No matter which light path is selected, the lever must be moved to the end (Figure 7).

(Figure 8)

4-10 Assemble photographic tube and photographic eyepiece

(1) Unscrew the dust cover and screw in the camera connector department.


(2) Insert the photographic eyepiece into the tee eyepiece holder.

(3) Put the camera tube on the tee eyepiece holder; then tighten the locking screw (Figure 8).

4-11 Adjusting the brightness of the transmitted illumination of the base

(1) Rotate the brightness adjustment knob as

shown by the icon ① on the base Button, turn clockwise, the brightness will increase; otherwise, the brightness will decrease (Figure 9)

(Figure 9)

5. Technical Parameters

5-1 Stereo microscope

Eyepiece	Standard configuration		Additional large objective lens					
			0.5X		1.5X		2X	
	Working distance100mm		Working distance165mm		Working distance45mm		Working distance30mm	
	Magnification	Field of view range	Magnification	Field of view range	Magnification	Field of view range	Magnification	Field of view range
10X/20	7X	28.6	3.5X	57.1	10.5X	19	14X	14.3
	45X	4.4	22.5X	8.9	67.5X	3	90X	2.2

5-2 Base electrical specifications

Mode	Composition
Power supply	220V-50Hz- 110V-50/60Hz
Transformer	Built-in 0-12V, continuous dimming
Lighting source	Ring 3W LED

★ This working distance is fixed and does not change with the change of magnification.

★ After using the auxiliary objective lens, total magnification = objective lens magnification × eyepiece magnification × auxiliary objective lens magnification

Eyepiece field of view

Object field of view (mm) = $\frac{\text{Eyepiece field of view}}{\text{Objective lens magnification} \times \text{auxiliary objective lens magnification}}$

★ The magnification of the film on the photographing device = the magnification of the objective lens (× the magnification of the auxiliary objective lens) × the magnification of the photographic eyepiece.

★ The magnification of the film on the camera device = the magnification of the objective lens (× the magnification of the auxiliary objective lens) × the intermediate magnification of the C-type camera tube.

6. Fault handling

If the performance of the microscope cannot be fully utilized due to unfamiliar use, the following table can provide you with some solutions.

6-1 General troubleshooting

Problem	Cause	Solution
1. Double images do not overlap	Incorrect interpupillary distance adjustment	Correct interpupillary distance
	Incorrect diopter adjustment	Re-adjust the diopter
	The left and right eyepieces have different magnifications	Install the same eyepiece
2. There is dirt in the field of view	There is dirt on the specimen	Clean specimen
	Dirt on the surface of the eyepiece	Clean the eyepiece
3. Unclear image	There is dirt on the surface of the objective lens	Clean the objective lens
4. Image is not clear when zooming	Incorrect diopter adjustment	Re-adjust the diopter
	Incorrect focus	Refocus
5. Focusing hand wheel is not flexible	The focus handwheel is locked too tightly	Loosen appropriately
6. During the observation process, the microscope body descends by itself, making the image unclear	Focusing hand wheel is too loose	After the focusing wheel is unscrewed, there are two spacers, which need to be close to the body, and then rotate the focusing wheel to tighten
7. There are cuts in the field of view of the observation eyepiece or photographic camera	The lever has not moved in place	Move the lever into place
8. When adjusting the focus, the image on the video monitor is not clear	Improper adjustment of camera focus	Use the focus ring on the camera tube to re-adjust the focus depth of the camera
9. Eyes tired easily	The diopter is not adjusted correctly	Correctly adjust the diopter
	Inappropriate lighting	Adjust the dimming knob
10. The bulb does not light up when the switch is turned on	No power	Check the connection of the power cord
	The bulb is not inserted	Insert correctly
	The light bulb is broken	Replace
11. The bulb burned out suddenly	Use of non-designated bulb	Replace with designated bulb
	Voltage is too high	Control voltage (if using voltage regulator)
12. Insufficient lighting	Use of non-designated bulb	Replace with designated bulb
	Voltage is too low	Increase voltage
13. The bulb is flickering or unstable	The light bulb is about to break	Replace
	Unstable lamp installation	Check and insert firmly

6-2 Video troubleshooting

Problem	Cause	Solution
1. The image on the display is cut	The lever has not moved in place	Move the lever into place
2. The image on the display has dirty spots	There is dirt on the specimen	Clean specimen
	There is dirt on the surface of the objective lens	Clean the objective lens
3. Image is not clear when zooming	The image is not clear when in a high magnification	Set the highest magnification, refocus

**ATTENTION!**

Before using this device, read this guide which contains important operating instructions for safe usage and control for compliance with applicable standards and regulations.

FCC Requirements:

•Products authorized under Part 15 using SDoC or Certification require a label containing one of the following compliance statements
(1)Receivers associated with licensed device service operations:

This device complies with part 15 of the FCC Rules. Operation is subject to the condition that this device does not cause harmful interference.

(2)Stand-alone cable input selector switch:

This device complies with part 15 of the FCC Rules for use with cable television service.

(3)All other devices:

•This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and

(2) this device must accept any interference received, including interference that may cause undesired operation.

CE Requirements:

•(Simple EU declaration of conformity) Hong Kong Svbon Technology Co.,Ltd declares that the equipment type is in compliance with the essential requirements and other relevant provisions of RED Directive 2014/30/EU and the ROHS Directive 2011/65/EU and the WEEE Directive 2012/19/EU; the full text of the EU declaration of conformity is available at the following internet address: www.svbon.com.

•Disposal

The crossed-out wheeled-bin symbol on your product, literature, or packaging reminds you that in the European Union, all electrical and electronic products, batteries, and accumulators (rechargeable batteries) must be taken to designated collection locations at the end of their working life. Do not dispose of these products as unsorted municipal waste. Dispose of them according to the laws in your area.

**IC Requirements:**

CAN ICES-3(B)/NMB-3(B)

Avoid Choking Hazard

Small Parts. Not for children under 3 years.

Approved Accessories

WARNING

•This device meets the regulatory standards when used with the Svbon accessories supplied or designated for the product.

•For a list of Svbon-approved accessories for your item, visit the following website: <http://www.Svbon.com>

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Warranty Card

Product Model	
Purchasing Date	
Defect Reason	
Dealer Name	
Telephone	
User's Name	
User's Address	
User's Email	

Remarks:

1. This guarantee card should be kept by the user, no replacement if lost.
2. Most new products carry a one-year manufacturer's warranty from the date of purchase.
3. The user can get warranty and after-sales service as below:
 - Contact the seller where you buy.
4. For warranty service, you will need to provide a receipt proof of purchase from the actual seller for verification

Exclusions from Warranty Coverage:

1. To any product damaged by accident.
2. In the event of misuse or abuse of the product or as a result of unauthorized alterations or repairs.
3. If the serial number has been altered, defaced, or removed.

MADE IN CHINA

Hong Kong Svbonny Technology Co.,Ltd

Add: Unit B, 5th Floor, Gallo Commercial Building, 114-118
Lockhart Road, Wanchai, Hong Kong

Facebook: [facebook.com/svbonny](https://www.facebook.com/svbonny)

E-mail: info@svbonny.com

Web: www.svbonny.com

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